



Kluba Medical

Medibino Neo

by Kluba Medical



Head positioning system for premature
and newborn babies for clinic use



Medibino Neo - Head Positioning System

- Medibino Neo is a head positioning system for premature and newborn babies for use in neonatology and obstetrics.
- The system consists of 3 components: gel ring, inner cushion and cotton cover
- Additionally, disposable covers are offered
- CE Medical Product Class I according to Directive 93/42/EEC December 2020

Application / Intended Use

- Stabilisation of the middle position of the head for premature babies < 33 weeks
- Stabilisation of head positioning during ventilation
- Pressure-relieving head support to prevent cranial deformations

Components of the Medibino Neo Head Positioning System



Soft gel ring

- Available in sizes S, M, L
- Size L: 130 mm x 110 mm x 25 mm
- Size M: 110 mm x 95 mm x 22 mm
- Size S: 90 mm x 80 mm x 18 mm
- Non-sterile, wipe disinfectable

Inner cushion made of 2.5-layer laminate, SMPTX 118g/sq m

- Available in sizes S, M, L to match the gel ring
- Waterproof and wipe disinfectable
- Filled with soft gel balls

Cotton cover

- Available in sizes S, M, L to match the gel ring
- Hotel closure with insert for ring and inner cushion
- Washable at 60 °C

All components are also available separately

Article Overview & Pricing (1/3)

Pos.	Item no.	Designation	Quantity / year	PCS	Price / piece
1	24021703-L	Medibino Neo gel ring size L	1,000 10,000 20,000 50,000	Piece	€ 35,00 € 26,00 € 23,00 € 20,00
2	24021703-M	Medibino Neo gel ring size M	1,000 10,000 20,000 50,000	Piece	€ 34,50 € 25,50 € 22,00 € 19,00
3	24021703-S	Medibino Neo gel ring size S	1,000 10,000 20,000 50,000	Piece	€ 34,00 € 25,00 € 21,50 € 18,50

Article Overview & Pricing (2/3)

Pos.	Item no.	Designation	Quantity / year	PCS	Price / piece
4	25021703-L	Medibino Neo inner cushion Size L	1,000	Piece	€ 8,00
			10,000		€ 7,00
			20,000		€ 6,00
			50,000		€ 5,50
5	25021703-M	Medibino Neo inner cushion Size M	1,000	Piece	€ 7,00
			10,000		€ 6,00
			20,000		€ 5,00
			50,000		€ 4,50
6	25021703-S	Medibino Neo inner cushion Size S	1,000	Piece	€ 7,00
			10,000		€ 6,00
			20,000		€ 5,00
			50,000		€ 4,50

Article Overview & Pricing (3/3)

Pos.	Item no.	Designation	Quantity / year	PCS	Price / piece
7	26021703-L	Medibino Neo Cover cotton Size L	1,000 10,000 20,000 50,000	Piece	€ 7,00 € 6,00 € 5,00 € 4,00
8	25021703-M	Medibino Neo Cover cotton Size M	1,000 10,000 20,000 50,000	Piece	€ 7,00 € 6,00 € 5,00 € 4,00
9	26021703-S	Medibino Neo Cover cotton Size S	1,000 10,000 20,000 50,000	Piece	€ 7,00 € 6,00 € 5,00 € 4,00



Kluba Medical GmbH
+49 176 219 70 684
n.klingen@kluba-
medical.com

*For what's worth
protecting the most
in life*